



GET IN THE
HEARING
LOOP

hearingloss.org

For Audiologists & Hearing Aid Providers

Please advise me if my hearing instrument(s) have telecoils. If they do, please activate and program them for use in a hearing loop. Also, please write any additional instructions for their use on the back of this card. If I don't have a telecoil(s), can it be retrofitted or is there an accessory remote control or streamer option available?

Program Recommendations

For open RIC fittings: provide a manually accessible T only program

For closed/occluded fittings: provide an M+T program for listening in social settings where it is important to hear my own voice, though lower Mic gain may be beneficial

For closed/occluded fittings: provide a manually T only program for listening in high noise environments like airports or train stations

Remember to verify Mic and telecoil responses to be sure they **match** in gain, output and frequency response in user programmed mode (use ANSI-SPLIV test procedure)

Programming questions?

Contact: Juliëtte Sterkens, AuD, jsterkens@hearingloss.org

I would like to use my telecoil(s) in a hearing loop or when a neckloop is available.

Please write clear instructions on how to use my manually accessible telecoil program in my hearing instruments.

For example: The numbered position of the program, number of beeps or voice prompt, or location of the switch to activate the t-coil.



Audiologist/hearing care provider name & contact information